



COVIDIEN

positive results for life™

Sharps Disposal Foot Pedal Carts

- Foot pedal provides easy one handed disposal
- Four wheels provide mobility
- Locking casters prevent movement during use
- Telescopic handle for convenience
- Foot pedal guard provides temporary lock in high traffic areas
- Universal key and lock provide security to prevent unwanted removal of container
- Temporary open button keeps lid open during procedures
- Weighted base keeps container upright throughout use
- Easy removal minimizes lifting



Mobile. Stable. Secure.

Slide Lid Foot Pedal Carts



8940FP (Container not included)



8938FP



8981FP

Hinged Lid Foot Pedal Carts



8941FP (Container not included)



8991FP



8980FP

Order Information

Foot Pedal Carts, Slide Lid

Code	Description	Ship Case
8940FP	SharpsCart™ molded cart for 12 & 18 gal containers	1
8938FP	SharpsCart™ wire cart for 12 & 18 gal containers	1
8981FP	SharpsCart™ wire cart for 8 gal containers	1

Foot Pedal Carts, Hinged Lid

Code	Description	Ship Case
8941FP	SharpsCart™ molded cart for 12 & 18 gal containers	1
8991FP	SharpsCart™ wire cart for 12 & 18 gal containers	1
8980FP	SharpsCart™ wire cart for 8 gal containers	1

OSHA Bloodborne Pathogen Standard (29 CFR 1910.1030 (d) (4) (iii) (A) (2)):

During use... containers for contaminated sharps shall be maintained upright throughout use.

NIOSH/CDC Guidelines. *Selecting, Evaluating, and Using Sharps Disposal Containers* (DHHS (NIOSH Publication No. 97-111))

P.2- Containers (including those designed to be kicked or wheeled), should be stable when placed on a horizontal surface and when used as described in the product labeling.



COVIDIEN, COVIDIEN with Logo and ™ marked brands are trademarks of Covidien AG or an affiliate. ©2008 Covidien AG or an affiliate. All rights reserved.

H6289 10M 08/08