

Safety Lancets









- Easy to use, pressure activated safety lancets
- Designed for both patients and healthcare professionals
- Needle retracts after use to leave no possibility of reuse
- New easy to read, bi-lingual packaging

Dual use 18G blade lancet is ideal for tests requiring larger blood samples. May be used on finger or heel.

High flow 21G lancet is used when a larger blood sample is needed. Used for testing blood gas, cholesterol and hemoglobin.

Medium flow 23G lancet is perfect for use when testing blood glucose.

Low flow 30G offers minimal pain as the thinnest lancet in the line.

18 Gauge <i>Dual Use</i> High Flow- Finger stick Low Flow- Heel stick	21 Gauge High Flow	23 Gauge Medium Flow	30 Gauge Low Flow
1.8 mm Depth Setting	2.2 mm Depth Setting	2.2 mm Depth Setting	1.6 mm Depth Setting
			
18G x 1.8mm Item # 28-1810 NDC # 86227-0181-05	21G x 2.2mm Item # 28-2110 NDC # 86227-0211-05	23G x 2.2mm Item # 28-2310 NDC # 86227-0231-05	30G x 1.6mm Item # 28-3010 NDC # 86227-0301-15
			

**For Ordering Information
Contact CHARTECH INC.
1-877-658-1849
Fax 770-442-2843
sales@chartech.net**